

Flexible range of medical-graded carts for clinical and research purposes



- ▶ Wide range of carts for any operator needs
- ▶ Compact and ergonomic design
- ▶ Compatible with most COSMED devices
- ▶ Electrified and non-electrified configurations
- ▶ Multiple supports, including shelves, PC, monitor and cylinder holders
- ▶ Quick and effective cleaning

COSMED features a wide range of medical carts to ensure that any configuration needs are met.

The ergonomic design, the high sliding wheels with locking system and the multiple supports and add-ons are conceived to aid the operator in the daily routine, from quick and effective cleaning to an easy transportation of devices and accessories.


Carts with built-in isolation transformer are provided with an electrified supporting column, representing the best solution to ensure high performance with no loss of functionality. Depending on the model, different shelves, cylinder holders, drawer and supports are included.

The latest X3 and X3+ generation represents the best choice of configurable medical carts whenever floor space is an issue, being compact, ergonomic and adjustable. Whether they are configured with an isolation transformer or in non-electrified setups, X3 and X3+ can be customized with many add-ons to meet any operator need.

All carts are compatible with most COSMED devices, including Quark Series, K5, Q-NRG and BOD POD GS-X.



1-cylinder holder (left) and 3-cylinder holder (right) carts



Basic configuration of X3 (left) and X3+ (right) carts

Technical Specifications

Cart	3-cylinder holder	1-cylinder holder	X3	X3+
REF	C02900-01-04 (230V)	C03550-01-04 (230V)	A-258-900-004	A-258-900-003
	C02900-02-04 (120V)	C03550-02-04 (120V)		
Dimensions (cm)	70x67x101	53x68x107	54x58x147	54x58x152
Max. weight (kg)	113	89	83	83
Max. load (kg)	47	47	46	46
Height adjustment	-	-	-	85-110 cm
Electrified	Yes	Yes	No	No
Power supply	220-240/100-120 Vac, 50/60 Hz		-	-
Packaging (LxWxH)	80x120x98 cm	80x120x95 cm	90x40x60 cm	78x80x110 cm
Compatible with	Quark Series, Quark C12x/T12x		Quark Series (GFC only), BOD POD GS-X, Q-NRG, K5, Quark C12x/T12x	
Options and accessories				
Isolation transformer	●	●	○	○
Articulated arm	○	○	○	○
Quark holder	●	●	○	○
Cylinder holder	● (3)	● (1)	○ (1 or 2)	○ (1 or 2)
Syringe holder	-	-	○	○
Double monitor VESA holder	○	○	○	○
PC holder	○	○	○	○
Shelf	●	●	○	○
Drawer	●	-	-	-
Basket	-	-	●	●
Standards				
	IEC 60601-1 (safety) IEC 60601-1-2 (EMC) MDR 2017/745	IEC 60601-1 (safety) IEC 60601-1-2 (EMC) MDR 2017/745	IEC 60601-1 (safety) IEC 60601-1-2 (EMC) MDR 2017/745	IEC 60601-1 (safety) IEC 60601-1-2 (EMC) MDR 2017/745

● Standard ○ Optional



COSMED Srl

Via dei Piani di Monte Savello 37
Albano Laziale - Rome 00041, Italy

+39 (06) 931-5492 Phone

+39 (06) 931-4580 Fax

info@cosmed.com | cosmed.com

Distributed by

